

# Pulp Disintegrator



## INTRODUCTION

Disintegration of a pulp sample is the mechanical treatment, in water, of interlaced fibers to separate them from one another.

## APPLICATIONS

Pulp

## SPECIFICATIONS

- Meets TAPPI T-205, CPPA C.4 and ISO 5263

## FEATURES

- Digital counter
- Includes lightweight (1.1 kg) polymer container
- 316 Stainless cover
- Integrated sensors for container proximity and cover lock enhance operator safety

## OPTIONS

- Hot disintegrator-heated stock container maintains temperature of pre-heated pulp samples

## INSTRUMENT SIZE

Depth: 521 mm (20.5 in.)  
Height: 425 mm (16.75 in.)  
Width: 305 mm (12 in.)

## WEIGHT

40 kg (87 lb)

## CONNECTIONS

Electrical: Specify voltage requirements when ordering

## ORDERING INFORMATION

Pulp Disintegrator	Catalog number 73-18
Hot Disintegrator (w/o container)	Catalog number 73-19
3 liter heated container	Catalog number 73-19-01
7 liter container	Catalog number 73-19-02

**The TMI Group  
of Companies**



[www.testingmachines.com](http://www.testingmachines.com)

**Testing Machines Inc.** 40 McCullough Drive, New Castle, DE 19720 USA  
Tel: 800-678-3221 (USA) or (302) 613-5600, Fax: (302) 613-5619, e-mail: [info@testingmachines.com](mailto:info@testingmachines.com)

**Büchel BV t/a Messmer Büchel** Fokkerstraat 24, 3905 KV Veenendaal, Netherlands  
Tel: +31 (0)318 521500, Fax: +31 (0)318 540358, e-mail: [buchel@testingmachines.com](mailto:buchel@testingmachines.com)

**Lako Tool** 7400 Ponderosa Road, Perrysburg, OH 43552 USA  
Tel: (419) 662-5256, Fax: (419) 662-8225

**Lawson Hemphill** 1658 G A R Highway, Swansea, MA 02777 USA  
Tel: (508) 679-5364, Fax: (508) 679-5396, e-mail: [information@lawsonhemphill.com](mailto:information@lawsonhemphill.com)

**Adamel Lhomargy SAS** Z.A. de l'Habitat, Batiment 6, Route d'Ozoir, 77680, Roissy en Brie, France  
Tel: +33 (0)1 64409210, Fax: +33 (0)1 64409211, e-mail: [adamellhomargy@testingmachines.com](mailto:adamellhomargy@testingmachines.com)

**TMI Canada** 40 McCullough Drive, New Castle, DE 19720 USA  
Tel: (302) 613-5600, Fax: (302) 613-5619 e-mail: [canada@testingmachines.com](mailto:canada@testingmachines.com)